



# An Update on the Timepix2 (It's Finely Here!)



# The Principal Acknowledgments

S. George<sup>1,2</sup>, P. Petty<sup>1,2</sup>, T. Campbell-Ricketts<sup>1,2</sup>,  
S. Wheeler<sup>3</sup>, D. Turecek<sup>1,6</sup>, L. Tlustos<sup>1,3,4</sup>, C. McLeod<sup>2,3</sup>

The CERN Medipix2 Design Team

(M. Campbell<sup>4</sup>, W. Wong<sup>4,5</sup>, X. Llopart-Cudie<sup>4</sup>,  
R. Ballabriga-Sune<sup>4</sup>, T. Poikela<sup>4</sup>, and J. Alozy<sup>4</sup>),  
E. Semones<sup>2</sup> & the NASA Space Radiation Analysis Group<sup>2</sup>  
And at HIMAC S. Kodaira<sup>8</sup> & H. Kitamura<sup>8</sup>

<sup>1</sup> Physics Department, University of Houston, Houston, TX, USA <pinsky@uh.edu>

<sup>2</sup> NASA SRAG, Johnson Space Center, Houston, TX, USA

<sup>3</sup> NASA, Johnson Space Center, Houston, TX, USA

<sup>4</sup> Medipix Project, CERN, Geneva, Switzerland

<sup>5</sup> XIE, Friburg, Germany

<sup>6</sup> Currently: University of Geneva, Geneva, Switzerland

<sup>7</sup> ADVACAM s.r.o., Prague, Czech Republic

<sup>8</sup> HIMAC, QST/NIRS, Inage (Chiba) Japan.

# Some Medipix Collaborations Genealogy

- ◆ **Medipix (1)**—Formed in the early 1990's and ended with the formation of Medipix2...

Medipix Chip—64 x 64 170  $\mu$ m pixel Photon Counting (PC)

- ◆ **Medipix2**—Formed in the late 1990's (**still active**)

Medipix2 MXR—256 x256 55  $\mu$ m 2-Threshold PC Frame (250 nm IBM)

**Timepix** (2006)—256 x256 55  $\mu$ m **TOT** or TOA & PC Frame (250 nm IBM)

**Timepix2** (2018)—256 x256 55  $\mu$ m **TOT+TOA** & PC Frame (130 nm TSMC)

- ◆ **Medipix3**—Formed in 2006 (**still active**)

Medipix3, Medipix3.1...


Medipix3RX (~2012)—256 x256 55  $\mu$ m pixel, Charge Summing PC

**Timepix3** (2014)—256 x256 55  $\mu$ m pixel **TOT+TOA** & PC Frame & **Data-Driven** Readout... (135 nm IBM)

- ◆ **Medipix4**—Formed in late 2016 (**just starting up**)


Medipix4 (PC) and **Timepix4** (TOT+TOA & PC) Frame & Data Driven (65 nm TSMC) 55  $\mu$ m pixel, but number TBD...

# Both the Timepix and Timepix2 are Hybrid Pixel Detectors Using the same Sensors



**Detector and electronics readout are optimized separately**  
**Bonding is done with the FlipChip® process...**

# Schematic Cross Section of a Hybrid Pixel Detector for Penetrating Charged Particles



# Common Digitization Methods

## Time-Over-Threshold (TOT)

This is a “Wilkison” type Analog to Digital Converter (**ADC**)

A “Pulse-Shaping” circuit regularizes the input current into a pulse whose pulse-height is proportional to the total charge collected over the shaping time.

The input voltage pulse is matched to a Threshold discriminator to output a digital “1” level for the duration of the time that the input is above the Threshold.

The input voltage is discharged by a constant current during which a clock is counted to determine the digital value for the input charge. (Note: the Timepix TOT counts any time the discriminator output is “1” and the frame is open...)

## Time-Of-Arrival (TOA)

The discriminator output described in the TOT method above, triggers a counter to begin counting as soon as the digital level changes to “1”.

The end of the Frame stops the TOA counters. (A Common Stop Time to Digital Converter (**TDC**)).

**The 100 MHz clock frequency limits the resolution in both Timepix and Timepix2 cases...**

# How are the Timepix and Timepix2 SIMILAR

## ◆ They have the SAME:

**Pixel footprint (55  $\mu$ m square).**

**256 x 256 square matrix.**

**Sensor chips and bump-bonding.**

**Can be connected to interfaces by wire-bonds (Although some pad placements differ.).**

**Both use a “Frame-Based” DAQ.**




# So, what does the Timepix2 ADD?

## ◆ The Timepix2 has the following additions:

- TSMC 130 nm CMOS Technology (Timepix uses GF/IBM 250 nm)
- Can be connected to interfaces with “Through Vias.”
- A 28-bit pseudo-random output register (compared to 14-bits)
- Simultaneous TOT **and** TOA mode (compared to TOT **or** TOA)
- Employs “Ghost” track rejection.
- Employs and increased optional Adaptive Gain frontend preamp
- Frontend designed to collect “holes.”
- Has 8 externally accessible “Digital-Pixels.”
- Direct access to Bias-Voltage “Backside Pulse.”
- Uses separate readout and digitization clocks






# Example Data

- ◆ Plots show example and ‘average’ clusters.
- ◆ Cut on cluster size for primary peak
- ◆ Average/ensemble clusters formed by stacking clusters on top of each other aligned by the center of mass (assumption - all clusters within cut are approximately similar’)
- ◆ Profile plots show ‘lines’ through average clusters (i.e. not a projection, where columns/rows are summed)
- ◆ Timepix2 tested with Adaptive Gain (AG) turned on and off.


# Ag mode max TOT ~2x non Ag mode

Si 350 MeV/A @ 30 Degrees, Ag on vs Off


Profile through cluster long axis




Profile through cluster short axis




Ag Off Ensemble Track




Ag On Ensemble Track




Ag Off Example Track



Ag On Example Track




Cluster size distribution



In this case 1.92x (undershoots simulation predictions, generally observed throughout data)


# Previous Simulation Results

Note TOT at saturation with AG on is approximately 2x saturation with AG Off




## Fe 290 MeV/A @ 0 Degrees, Comparison to Timepix


Profile through cluster long axis




Profile through cluster short axis




Timepix1 Ensemble Track




Timepix2AG Ensemble Track




Timepix1 Example Track



Timepix2AG Example Track




Cluster size distribution




# Ag mode effects cluster skirt shapes

Si 350 MeV/A @ 75 Degrees, Ag on vs Off


Profile through cluster long axis




Profile through cluster short axis




Ag Off Ensemble Track




Ag On Ensemble Track




Ag Off Example Track



Ag On Example Track




Cluster size distribution




**Skirt larger and different shape with Ag off**

# Fe 290 MeV/A @ 80 Degrees, Comparison to Timepix


Profile through cluster long axis




Profile through cluster short axis




Timepix1 Ensemble Track




Timepix2AG Ensemble Track




Cluster size distribution



Timepix1 Example Track



Timepix2AG Example Track




Here when compared to reference Timepix1 skirt is larger




# Here when compared to reference Timepix1 skirt is the same size (much lower LET)

N 430 MeV/A @ 75 Degrees, Comparison to Timepix


Profile through cluster long axis




Profile through cluster short axis




Timepix1 Ensemble Track




Timepix2AG Ensemble Track




Timepix1 Example Track



Timepix2AG Example Track



Cluster size distribution




# Comparison of front end to Timepix1


- ◆ Tested Timepix and Timepix2 both had 500um sensors at 100V.
- ◆ In principal average clusters should be directly comparable.
- ◆ Can test Timepix2 by comparing to the ADC counts to the calibrated Timepix energy at the same cluster position...

# N 430 MeV/A @ 30 Degrees, Comparison to Timepix


Profile through cluster long axis




Profile through cluster short axis




Timepix1 Ensemble Track




Timepix2AG Ensemble Track




Timepix1 Example Track



Timepix2AG Example Track



Cluster size distribution





Photon calibration:




where  $t_{\text{thr}}$  is time over threshold,


and  $E$  is energy,  
becomes inaccurate above  $\sim 100$  keV.

Pixel behavior is still unclear above  
 $\sim 1400$  keV.

Cluster Volume 5 MeV Hydrogen at 60° polar angle




Cluster Volume 5 MeV Hydrogen at 60° polar angle



# Ghost Tracks

- ◆ "Ghost tracks" occur when the shutter overlaps a high TOT signal. This either cuts off the TOT signal prematurely, or counts ongoing TOT inputs due to events occurring just prior to the Frame opening
- ◆ Gives the cluster a 'decapitated' appearance, and results in incorrect measurement of TOT.
- ◆ For ISS, with lots of stopping protons effectively limits the minimum frame time to 10 mS, and the rate in an isotropic field to a few kHz depending on how many overlaps are acceptable.
- ◆ Timepix2 has special logic to stop this problem occurring. In principal this means we can now exploit one of the principal advantages of pixels - their high rate capabilities!





# Timepix2 Olympics @ HIMAC

- ◆ Previous testing with a Neon beam measuring very short frame lengths has indicated that the 'ghost track' suppression logic works. At the most recent HIMAC run we decided to exercise the Timepix2 under high particle fluxes.
- ◆ Beam starting at 100 particles per spill (PPS). Measured up to  $1\text{E}8$  PPS. Timepix2 measurement does not exactly correspond to this, as some particles miss Timepix2.
- ◆ Frame length set to measure single clusters (0.1 s at 100 PPS) and then scaled down proportional to the PPS.
- ◆ So at  $1\text{E}8$  PPS we had a 100 nS frame length (it was not clear if the Advadac readout system would handle frame lengths this short correctly, but no problems were observed)




# Results




Error bar shows beam FWHM of beam

# Ghost-Track Suppression

Maximum Cluster TOT  
(Cluster height)




Total Cluster ADC Counts



350 MeV/A Si 28 beam

**Timepix2 performs well for very short frame lengths**

# Backside Pulse Observed With a Scope at HIMAC



Backside pulse could probably be a bit more sensitive, signal was barely above noise for 100 MeV Ne strike at 0 degrees incident angle. Desired range is about 7 M electrons to about 300 M electron (higher would be nice!) S. George

## Continuing Efforts

- ◆ Harnessing the *Backside Pulse* using one of more of the Digital Pixels to allow reasonable estimate of very high LET events
- ◆ Devising an economical *Calibration* routine for the Adaptive-Gain frontend.
- ◆ Settling on an optimal set of Operational *DAC values* and *supply voltage* (VDD & VDDA) for the 130 nm CMOS frontend.




**Questions...?**

ADC counts (not calibrated)

Timepix 2 ADC Counts

Timepix 1 vs 2 Calib Data



Inflection point may be due to Timepix 'advanced' calibration which starts at 950 keV for this assembly


Limit of Timepix 1 calibration/start of volcano (2.2 MeV)

**Different markers show different polar angles**


**Tpx2 appears to measure monotonically up to at least 2 MeV**



# Voltage-dependent feedback capacitance







$$\frac{1}{\dots}$$



When the transistor is active,  
the feedback capacitance is  
~125fF  
(Dimensions transistor 1/10)


**M. Manghisoni et al. "Dynamic Compression of the Signal in a Charge Sensitive Amplifier: from Concept to Design" IEEE Trans. Nucl. Sci. Vol. 62, No. 5, 2015**



Amplitude, ToT and ToA for the front end with adaptive gain ON, front end power supply 1.2V

(a low power mode for the discriminator is also considered and has an impact on the ToA (discriminator propagation delay))

(Gain: ~20mV/ke- (CBEST), ~16mV/ke- (CWORST))




# Pixel Digital Operating Modes


**28 bits/pixel:**  
**2 x 4b chains**  
**2 x 10b chains**

Mode Name	Description	1 <sup>st</sup> Counter	2 <sup>nd</sup> Counter
<b>ToT10/ToA18</b>	Simultaneous ToT and 1 <sup>st</sup> hit ToA*	10-bit ToT	18-bit ToA
<b>ToT14/ToA14</b>	Mode options (programmable): 1) 1 <sup>st</sup> hit or integral ToT 2) Overflow (wraparound) of ToA counter	14-bit ToT	14-bit ToA
<b>ContToT10/Event4</b>	Continuous read/write ToT  Mode options: 1) 1 <sup>st</sup> hit or integral ToT (programmable) 2) Supplementary 4-bit eventing counting (readout optional)	10-bit ToT 4-bit #Events	10-bit ToT 4-bit #Events
<b>ContToT14</b>	Continuous read/write ToT  Mode option: 1 <sup>st</sup> hit or integral ToT (programmable)	14-bit ToT	14-bit ToT
<b>ContToA10</b>	Continuous read/write 1 <sup>st</sup> hit ToA	10-bit ToA	10-bit ToA
<b>ContToA14</b>		14-bit ToA	14-bit ToA
<b>ContEvent10</b>	Continuous read/write event counting	10-bit #Events	10-bit #Events
<b>ContEvent14</b>		14-bit #Events	14-bit #Events


Mode Name	Description	1 <sup>st</sup> Counter	2 <sup>nd</sup> Counter
ToT10/ToA18	Simultaneous ToT and 1 <sup>st</sup> hit ToA*  Mode options (programmable): 1) 1 <sup>st</sup> hit or integral ToT 2) Overflow (wraparound) of ToA counter	10-bit ToT	18-bit ToA
ToT14/ToA14		14-bit ToT	14-bit ToA




Mode Name	Description	1 <sup>st</sup> Counter	2 <sup>nd</sup> Counter
ToT10/ToA18	Simultaneous ToT and 1 <sup>st</sup> hit ToA*	10-bit ToT	18-bit ToA
ToT14/ToA14	Mode options (programmable): 1) 1 <sup>st</sup> hit of integral ToT 2) Overflow (wraparound) of ToA counter	14-bit ToT	14-bit ToA



Mode Name	Description	1 <sup>st</sup> Counter	2 <sup>nd</sup> Counter
<b>ContToT10/Event4</b>	Continuous read/write ToT Mode options: 1) 1 <sup>st</sup> hit or integral ToT (programmable) 2) Supplementary 4-bit eventing counting (readout optional)	10-bit ToT 4-bit #Events	10-bit ToT 4-bit #Events
<b>ContToT14</b>	Continuous read/write ToT Mode option: 1 <sup>st</sup> hit or integral ToT (programmable)	14-bit ToT	14-bit ToT





# Power Consumption in the Matrix

- **Analogue:**
  - 5  $\mu\text{A}$ /pixel @ 1.2 V
  - Unused pixels can be powered down to a few nA
  - With all 65536 pixels powered up: <400 mW
- **Digital:**
  - 0.24 mW/superpixel during open Shutter in simultaneous ToT/ToA mode with 100 MHz ToTClk and 100 MHz ToAClk, assuming sparse activity and superpixel clock gating enabled
  - ToT and ToA clocks can run at lower frequencies
    - ToAClk is output from a clock divider taking ToTClk as input
  - The clocks in a double column are disabled when all pixels in the double column are masked
  - Full matrix digital power with sparse activity and 100 MHz clocks: <500 mW
- **Total: <900 mW for the full matrix assuming sparse hits and both clocks @ 100 MHz**
  - ~2 W for the full matrix in cases of full occupancy (not relevant to dosimetry)






# Maximum Frame Rates

## With a 100 MHz DCLOCK\_IN:

[illegible]



# Timepix2 v. Timepix

1. Improved Front-End TOT range and stability...
2. Simultaneous TOT and TOA (10 ns)...
3. Stable performance with 100 MHz clock...
4. Stable and linear internal digital and analog test pulses...
5. TSV capable (no wirebonds & secure mounting)...
6. 28-bit (Fungible) output buffer...
7. Improved DAC stability...
8. Multiple continuous operational modes...
9. Can use Timepix sensors...
10. **Bonus—8 Digital Pixels to couple external devices**